

## SKINTACT Cool Contact • Electrosurgical Neutral Electrodes

“SKINTACT Cool Contact Electro-surgical Neutral Electrodes provide safety of the highest degree. They are compatible with all common generators on the market and fully comply with relevant clauses of EN 60601-2-2:2009 guaranteeing patient safety and your peace of mind.”
























# SKINTACT Cool Contact • Electrosurgical Neutral Electrodes

The professional's choice


- Unique C-LINE manufacturing produces ultra soft hydrogel which conforms perfectly to the body contours providing maximum safety
- Low skin contact impedance ensures the neutral electrode stays as cool as possible during use
- Flexible low profile foam, easy to mould provides excellent security
- An enhanced adhesive border provides extra security should the neutral electrode be pulled inadvertently or come into contact with prepping solutions
- Available in a comprehensive range of sizes, adult/pediatric, split/unsplit and with or without pre-attached cable

## SKINTACT Cool Contact Electrosurgical Neutral Electrodes



							
<b>Order Code</b>	<b>RS05</b>	<b>RS06</b>	<b>RS03</b>	<b>RS25</b>	<b>RS26</b>	<b>RS27</b>	<b>RS23</b>
<b>Type</b>	unsplit	unsplit	unsplit	split	split	split	split
<b>Size</b>	adult	pediatric	neonatal	adult	pediatric	adult	neonatal
<b>Conductive Area</b>	134 cm <sup>2</sup>	76 cm <sup>2</sup>	38 cm <sup>2</sup>	118 cm <sup>2</sup>	66 cm <sup>2</sup>	120 cm <sup>2</sup>	34 cm <sup>2</sup>
<b>Weight Limit</b>		5 - 15 kg (11 - 33 lbs)	0.45 - 4.99 kg (1 - 11 lbs)		5 - 15 kg (11 - 33 lbs)		0.45 - 4.99 kg (1 - 11 lbs)
<b>Power Setting</b>		< 120 W	< 0.35 Amp		< 120 Watts		< 0.35 Amp
<b>Packaging Standard</b>	1 or 5 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 or 5 / pouch 50 / box	1 / pouch 50 / box	1 or 5 / pouch 50 / box	1 / pouch 50 / box
<b>with 3 m pre-attached disposable cable*</b>							
<b>Order Code</b>	<b>RS05B30</b>	<b>RS06B30</b>	<b>RS03B30</b>	<b>RS25A30</b>	<b>RS26A30</b>	<b>RS27A30</b>	<b>RS23A30</b>
<b>Connector</b>	Non-REM	Non-REM	Non-REM	REM	REM	REM	REM
<b>with 3 m pre-attached disposable cable*</b>							
<b>Order Code</b>	<b>RS05D30</b>	<b>RS06D30</b>	<b>RS03D30</b>	<b>RS25D30</b>	<b>RS26D30</b>	<b>RS27D30</b>	<b>RS23D30</b>
<b>Connector</b>	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack
<b>Packaging Standard</b>	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box

**Compatibility:** In any case, split Skintact Neutral Electrodes are compatible with the dynamic contact quality monitoring systems REM™, NESSY™ and ARM™.

## SKINTACT Cool Contact Reusable Cables\* and Adaptors

						
<b>Order Code</b>	<b>NEK300A</b>	<b>NEK300B</b>	<b>NEK300D</b>	<b>NEK300E</b>	<b>NEK300G</b>	<b>NEK300H</b>
<b>Connector (compatible)</b>	REM (Valleylab)	Non-REM (Valleylab)	6.3 mm jack (Standard)	4 mm banana	2 Pin adjustable (Eschmann)	6.3 mm jack (Erbe)
<b>Cable</b>	3 m	3 m	3 m	3 m	3 m	3 m

\* SKINTACT cables can also be supplied in 5 metre.

<b>Order Code</b>		<b>NEKA-AD</b>		<b>NEKA-AH</b>
<b>Description</b>		Adaptor A/B-D		Adaptor A-H

SKINTACT Cool Contact Electrosurgical Neutral Electrodes fully comply with relevant clauses of EN 60601-2-2:2009.

**SKINTACT Cool Contact Reusable Dispersive Cable with stamped Lot Number for easy tracking (types REM and NON-REM).**

